

<ul style="list-style-type: none"> <li>● egnet</li> <li>✗ ikke egnet</li> <li>▼ begrænset resistens</li> <li>○ ingen oplysninger</li> </ul>	Membrantype										
	Celluloseacetat	Mixed Cellulose Ester	Nylon	Polycarbonat	Polyethersulfon	Polypropylen	PTFE	PTFE, lamineret	PTFE, hydrofil	PVDF	Regeneret Cellulose
Tabellen gælder for løse membranfiltre. Er membranen indbygget i et filterhus, skal dettes resistens også tages i betragtning.											
Acetic Acid 1%	●	●	●	●	●	●	●	●	●	●	●
Acetic Acid 5%	●	●	●	●	●	●	●	●	●	●	●
Acetic Acid 10%	○	●	●	✗	●	●	●	●	●	●	●
Acetic Acid, Glacial	✗	✗	✗	✗	●	●	●	●	●	●	●
Acetone	✗	✗	●	▼	✗	●	●	●	●	✗	●
Acetonitrile	✗	✗	●	▼	●	●	●	●	●	●	●
Ammonium Hydroxide, 6N	✗	●	✗	✗	●	●	●	●	●	✗	✗
Amyl Alcohol	●	✗	●	●	✗	●	●	●	●	●	●
Amyl Acetate	○	✗	●	●	●	●	●	●	●	✗	●
Aniline	✗	✗	●	✗	▼	○	●	●	●	○	●
Benzene	●	●	●	▼	●	▼	●	●	●	●	●
Benzyl Alcohol	▼	●	●	●	✗	●	●	●	●	●	●
Boric Acid	●	●	▼	●	●	●	●	●	●	●	○
Boric Acid 5%	●	●	▼	●	●	●	●	●	●	●	○
Bromoform	✗	●	●	●	✗	●	●	●	●	○	○
Butyl Alcohol	●	●	●	●	●	●	●	●	●	●	○
Butyl Acetate	✗	✗	✗	●	●	●	●	●	●	○	●
Butyl Cellosolve	▼	✗	●	▼	●	●	●	●	●	○	○
Carbon Tetrachloride	▼	●	●	●	●	▼	●	▼	●	●	●
Cellosolve	●	✗	●	●	✗	●	●	●	●	●	●
Chloroform	✗	●	●	●	✗	▼	●	▼	●	●	●
Cottonseed Oil	●	●	●	●	●	●	●	●	●	○	○
Cyclohexane	●	●	●	●	●	●	●	●	●	○	●
Cyclohexanone	✗	✗	●	▼	✗	●	●	●	●	✗	●
Diethyl Acetamide	✗	✗	●	✗	✗	✗	●	●	●	○	●
Dimethyl Formamide	✗	✗	●	✗	✗	●	●	●	●	✗	▼
Dimethyl Sulfoxide (DMSO)	✗	✗	▼	✗	✗	●	●	●	●	✗	●
Dioxane	✗	✗	●	✗	✗	●	●	▼	●	▼	●



<ul style="list-style-type: none"> <li>● egnet</li> <li>✗ ikke egnet</li> <li>▼ begrænset resistens</li> <li>○ ingen oplysninger</li> </ul>	Membrantype										
	Celluloseacetat	Mixed Cellulose Ester	Nylon	Polycarbonat	Polyethersulfon	Polypropylen	PTFE	PTFE, lamineret	PTFE, hydrofil	PVDF	Regeneret Cellulose
Tabellen gælder for løse membranfiltre. Er membranen indbygget i et filterhus, skal dettes resistens også tages i betragtning.											
Ethyl Alcohol, >80%	●	▼	●	●	●	●	●	●	●	○	○
Ethyl Alcohol, <80%	●	●	●	●	●	●	●	●	●	●	●
Ethyl Ether	▼	▼	●	●	▼	●	●	●	●	●	●
Ethylene Dichloride	▼	▼	●	●	✗	●	●	●	●	○	○
Ethylene Glycol	●	▼	●	●	●	●	●	●	●	●	●
Formaldehyde	▼	✗	●	●	●	●	●	●	●	●	●
Formic Acid 50%	○	▼	○	●	○	○	●	●	●	○	○
Freon TF	●	●	●	●	●	●	●	●	●	●	○
Gasoline	●	●	●	●	●	●	●	●	●	●	●
Glycerine (Glycerol)	●	●	●	●	●	●	●	●	●	●	●
Hexane	●	●	●	●	●	●	●	●	●	●	●
Hydrochloric Acid, Conc.	✗	✗	✗	●	●	●	●	●	●	●	✗
Hydrochloric Acid, 6N	▼	✗	✗	●	●	●	●	●	●	●	○
Hydrochloric Acid, 12N	✗	✗	✗	●	●	●	●	●	●	●	○
Hydroflouric Acid	✗	✗	✗	●	●	●	●	●	●	●	○
Hydroflouric Acid, 10%	✗	✗	✗	●	○	●	●	●	●	○	○
Hydroflouric Acid, 35%	✗	✗	✗	●	○	○	●	○	●	○	○
Hydrogen Peroxide (30%)	●	●	●	●	●	●	●	●	●	●	○
Isobutyl Alcohol	●	●	●	●	●	●	●	●	●	●	○
Isopropyl Acetate	✗	✗	●	●	●	●	●	●	●	●	●
Isopropyl Alcohol	●	▼	●	●	○	○	●	●	●	●	○
Isopropyl Ether	●	●	○	●	○	○	●	●	●	○	○
Kerosene	●	●	●	●	●	●	●	●	●	●	○
Kodak KMER, FTFR	✗	✗	●	●	✗	●	●	●	●	○	○
Methanol	●	✗	●	●	●	●	●	●	●	●	●
Methyl Acetate	✗	✗	●	✗	✗	●	●	●	●	●	●
Methyl Ethyl Ketone (MEK)	✗	✗	●	▼	✗	●	●	●	●	○	●

<ul style="list-style-type: none"> <li>● egnet</li> <li>✗ ikke egnet</li> <li>▼ begrænset resistens</li> <li>○ ingen oplysninger</li> </ul>	Membrantype										
	Celluloseacetat	Mixed Cellulose Ester	Nylon	Polycarbonat	Polyethersulfon	Polypropylen	PTFE	PTFE, lamineret	PTFE, hydrofil	PVDF	Regeneret Cellulose
Tabellen gælder for løse membranfiltre. Er membranen indbygget i et filterhus, skal dettes resistens også tages i betragtning.											
Methyl Cellosolve	▼	▼	●	✗	●	●	●	●	●	○	○
Methyl Isobutyl Ketone	✗	✗	●	▼	✗	●	●	●	●	○	●
Methylene Chloride	✗	✗	▼	✗	✗	●	●	●	●	●	●
Nitric Acid, Conc.	✗	✗	✗	●	✗	✗	●	●	✗	●	✗
Nitric Acid, 6N	▼	●	✗	●	▼	▼	●	●	●	●	✗
Nitric Acid, 12N	✗	✗	✗	●	✗	○	●	●	▼	●	✗
Nitrobenzene	✗	✗	●	●	✗	●	●	●	●	●	●
Peanut Oil	●	●	●	●	●	●	●	●	●	○	○
Pentane	●	●	●	●	●	▼	●	●	●	●	●
Perchloric Acid 60%	✗	●	○	✗	○	○	●	●	●	○	▼
Perchloroethylene	●	●	●	●	✗	●	●	●	●	●	●
Phenol	✗	●	●	✗	○	○	●	●	●	○	○
Phosphoric Acid 85%	✗	●	○	✗	○	○	●	●	●	○	✗
Potassium Hydroxide, 6N	✗	✗	●	✗	●	●	●	●	●	○	✗
Propanol	●	●	●	●	●	●	●	●	●	●	○
Pyridine	✗	✗	●	●	✗	●	●	●	●	●	●
Shipley (AS-111, 340, 1350)	✗	✗	●	●	✗	●	●	●	●	○	○
Silicone	●	●	●	●	●	●	●	●	●	●	○
Sodium Hydroxide, 6N	✗	✗	✗	✗	●	●	●	●	●	○	✗
Sulfuric Acid, Conc	✗	✗	✗	✗	✗	✗	●	●	✗	✗	✗
Sulfuric Acid, 6N	▼	●	✗	●	▼	▼	●	●	●	●	✗
Tetrachloroethylene	●	●	○	●	○	○	●	●	●	○	○
Tetrahydrofuran	✗	✗	▼	✗	✗	▼	●	▼	●	●	●
Toluene	●	●	●	▼	✗	▼	●	●	●	●	●
Trichloroethane	▼	●	●	✗	✗	●	●	●	●	●	●
Trichloroethylene	●	●	●	✗	✗	●	●	▼	●	●	●
Triethylamine	●	●	●	▼	✗	●	●	●	●	○	○
Turpentine	●	●	●	●	●	●	●	●	●	○	○
Xylene	●	●	●	●	✗	▼	●	●	●	●	●